




Standard Operating Procedure for SEM3

[Note the locations of the various SEM3 components in this [Video: SEM3-00](#)]




*** Videos have no sound ***

Before beginning, ensure you have an active [reservation](#) for SEM3

1 - Loading a Sample:

- a) If *Sample Exchange Window* is not open, click  to open it [[Video: SEM3-1a](#)]
- Click *Chamberscope* tab
 - If  icon appears in the Chamberscope image, click to un-pause
 - Click  to initiate the vent cycle
 - Confirm by pressing the vent button in the pop-up window
 - A green status bar will appear at the bottom of the *Sample Exchange Window*... wait for the status bar to disappear, indicating completion of the vent cycle
- b) Gently slide open chamber door; place your sample stub in one of the holes in the stage. Make sure your sample is securely mounted to the stub (e.g. with vacuum compatible conductive tape)

Note: If loading multiple samples, limit stubs to one type per session. For example, do not mix flat stubs with 45° stubs or 90° stubs. This is to avoid collisions with the pole piece due to different sample heights. [[Video: SEM3-1b](#)]

- c) Watch the live CCD chamber view and gently slide the chamber door closed [[Video: SEM3-1c-e](#)]
-  Make certain the sample will not contact the pole piece (your sample should be below the 10mm mark on screen)
- d) Fully close the chamber door then click  in the *Sample Exchange Window*
- e) Wait for the pumping status bar to reach the 'full' mark and then disappear, indicating the pump-down sequence is complete (~ 5min). Additionally, the chamber/column vacuum status icon at the bottom of the main UI window will appear fully green: 


2 - Preparing to Image:

- a) At the top of the UI window, click **Stage** and select **Take Nav-Cam Photo**. Wait for the Nav-Cam operation to complete (~ 1min). An image of your sample should appear in Quadrant 3 (Q3, the lower-left quadrant). [[Video: SEM3-2a](#)]




- b) Lower magnification to its minimum by rotating the Magnification knob on the control panel counter-clockwise. [\[Video: SEM3-2b-c\]](#)
- c) Navigate to your sample by double-clicking a location in the Nav-Cam photo; the x-y stage will move to the specified location
- d) In the *Sample Exchange Window*, raise the sample stage to the 10mm marker by holding the middle mouse button and dragging upward; release the mouse button when your sample reaches the green 10mm marker [\[Video: SEM3-2d\]](#)




If your sample is magnetic, check the **Magnetic Sample** option







- e) Click Mode Selection:  and select **Mode 2: Optiplan** in the main UI [\[Video: SEM3-2e-i\]](#)
- f) Set accelerating voltage and emission current in the drop-down menus; recommended initial values are 2kV and 25pA:


NOTE: For those accustomed to using SEM1 in SMIF, you may be tempted to use high voltages (e.g. 20kV – 30kV), with the intent of getting high-resolution images. The Apreo does not require high voltages for high resolution imaging; 1kV – 5kV is typical.

- g) Click the Beam Control button on the side pane at the right of the UI: 
- h) Click  to turn on the beam; the button will turn from grey to orange and you will hear a pneumatic valve actuate within the SEM column
- i) Click  to zero the beam shift (this centers the beam within the column)

3 - Basic Imaging:

- a) Click within Quadrant 1 (Q1, upper left) to activate the quadrant; the data bar at the bottom of Q1 should be blue, indicating Q1 is active [\[Video: SEM3-3a-c\]](#)
- b) Click the Detectors button on the side pane 
 - Select **ETD**: (this assigns the Everhart-Thornley detector to the active quadrant, Q1)
 - Select **Secondary Electrons**:

- c) Click the green pause button in Q1 to turn on the detector and begin imaging: 
- Note: you can press the grey/black pause button at the top of the UI at any time to pause or un-pause image acquisition for the *active* quadrant: 
- d) Adjust Brightness and Contrast manually by rotating the knobs on the control panel -OR- perform the Auto Contrast & Brightness function by pressing:  [\[Video: SEM3-3d-f\]](#)
- e) Adjust Focus knob, then increase magnification to at least 5,000x and re-focus
- f) Click the **Link Z to FWD** icon: 
- With the sample in focus, pressing this icon will *link* the position of the sample to the focal length of the microscope... establishing the *Free Working Distance*, which is the distance between the sample and the bottom of the pole piece
 - Before **Link Z to FWD** is clicked, the icon appears with a question mark (as shown above); after clicking, the icon changes to: , indicating successful linking
 - After linking, the stage Z-coordinate will display the working distance between the pole piece and portion of the sample used during linking.
-  Attention: If your sample is not flat -OR- is mounted such that it is not parallel to the sample stage, you risk collision with the pole piece at smaller working distances. For *non-flat* samples, focus at the highest point on the sample and then link Z to FWD.
- g) You can now navigate across your sample to locate features of interest. Navigation can be accomplished in several ways: [\[Video: SEM3-3g\]](#)
- Mouse: Double-click a point on the Nav-Cam image OR the electron detector image to move the stage to that location. Note: Nav-Cam navigation is typically used at lower magnifications, to navigate to approximate locations on a sample (due to resolution limits of the optical Nav-Cam image)
 - Arrow keys: While imaging the sample with the ETD in Q1 (note: Q1 must be *active*), pressing the left / right / up / down arrow keys on the keyboard will shift the stage in each direction. The magnitude of the shift is 80% the field of view.

- Stage coordinates: Click the Stage button on the side pane: . You can enter in specific x- and y-coordinates to drive the stage to that location. Sample Rotation and Sample Tilt are also available



Attention: The use of Tilt should be done with caution. Tilting samples greatly increases the risk of a pole piece collision. Always view the LIVE CCD image when tilting. Tilt with caution! MAXIMUM tilt angle allowed is 30°

h) Some notes on Working Distance: [\[Video: SEM3-3h\]](#)

- At this point, your Working Distance (WD, as displayed in the data bar) should be *approximately* 10mm (you positioned your sample at the 10mm marker in the **Preparing to Image** section above)
- You should enable the live CCD view whenever you change WD (the z-coordinate of the sample) OR whenever you Tilt; this is to ensure there is NO chance of a pole piece collision
- A good WD for the ETD detector is 7mm
- A good WD for Immersion mode is 5mm
- The MINIMUM WD SMIF allows is 3.0mm; Do NOT use WD smaller than 3.0mm




Attention: If you have initiated stage movement and you fear a pole piece collision is imminent, you can halt stage movement by pressing the ESC key on the keyboard. Be ready to press ESC!

4 - Improving Image Quality:

a) Iteratively adjust Focus → X-Stigmation → Y-Stigmation to improve image [\[Video: SEM3-4a\]](#)

b) Perform Lens Alignment: [\[Video: SEM3-4b\]](#)


- Make note of some feature in the live image in Q1
- Click the Lens Alignment button: 
- The image will begin to wobble. Drag the crosshairs that appear on the wobbling image in the x- and y- directions. The goal is to minimize the wobbling in the image.
- After minimizing wobbling, click the lens alignment button again

5 - Capturing Images:

- a) Users can adjust frame averaging, dwell time and scan resolution by adjusting settings from the

dropdown menus at the top of the UI:  [Video: SEM3-5a]



- b) -OR- users can apply one of three preset scans available by pressing the **s1**, **s2**, **s3** buttons at the

top of the UI:  (in this image, S1 is selected) [Video: SEM3-5b]








- **s1** – TV scan rate with a 0.1s frame acquisition time (use for live imaging)
- **s2** – HDTV scan with a 0.2s frame acquisition time (use for reduced area live imaging)
- **s3** – this button can be edited / customized by the user; right click on the s3 button and select edit (use for image capturing by pressing **Ctrl + s3**)

Note: please do not edit s1 or s2... only edit s3

- c) There are two preset image capture buttons; pressing either of these will scan and capture the image and prompt the user to save the image file: [Video: SEM3-5c]

-  The **Snapshot** icon captures a ~1.5MP image in ~ 24 seconds
-  The **Photo** icon captures a ~6.5MP image in ~ 60 seconds




6 - Sample Removal / Ending Your Session:


- a) Click  and set Scan Rotation to 0° if not already [Video: SEM3-6a-f]
- b) Click  and set WD to 10mm or greater; If Tilt or Rotation were used, set both to zero; set stage translation X=0 and Y=0
- c) Click  at the top of the UI to zero beam shift
- d) Click  and then  to turn off the beam
- e) Click  and confirm by pressing vent in the pop-up window
- f) If Q3 was changed, set back to Nav-Cam
- g) If Nav-Cam zoom was used, click **View > Undo Digital Zoom** at the top menu of the UI
- h) If Q4 was changed, set back to Camera
- i) Wait for vent cycle to complete, then slide open chamber and retrieve your sample
- j) Slide chamber closed and press  [Video: SEM3-6i-j]
- k) End your [reservation](#) if necessary

- **APPENDIX 1 - Utilizing Multiple detectors (ETD, T1, T2, T3):** [[Video: SEM3-A1](#)]

- In addition to the ETD described above, the Apreo is equipped with a set of *Trinity* detectors, located within the column: T1, T2 and T3
- Following the procedure described in the **Basic Imaging** section above, users can assign different detectors to different quadrants (Q1, Q2, Q3, Q4)
- Users can opt between any of the detectors, depending on imaging needs
- Simultaneous imaging of up to 4 detectors (one in each quadrant) is possible


Note: Holding **Shift** and pressing a button will perform that action on all imaging quadrants:

- Shift +  captures a snapshot from all imaging quadrants
- Shift +  (un)pauses all imaging quadrants
- Shift +  performs auto brightness/contrast on all imaging quadrants




- **APPENDIX 2 - Immersion Mode:** [[Video: SEM3-A2](#)]

In the **Preparing to Image** section, **Optiplan** was selected as the imaging mode. Optiplan is suitable for a variety of samples and is capable of high resolution imaging. When ultra-high resolution is desired, users can enter **Immersion Mode** as follows:

- Verify your sample is NOT magnetic (if it is, you should have checked the **Magnetic Sample** in the **Preparing to Image** section above); do NOT proceed with immersion mode if your sample is magnetic!
- Ensure the CCD view is live (i.e. no green pause button in the upper left of Q4)
- Click the Stage button on the side pane: .
- Ensure that you have performed the Link Z to FWD operation (as described in the **Basic Imaging** section above); if linking was performed, there should be a down-pointing arrow beside the z-coordinate:



If not, perform Link Z to FWD as described in the **Basic Imaging** section

- With the CCD image live, change your working distance to 5.0mm (this is a good starting point for immersion mode)
- Click the down arrow on the Mode Selection icon:  and select **Mode 3: Immersion**
- Assign your detectors (T1, T2, T3) to your desired quadrants, as described in the **Basic Imaging** section; Note: The ETD is not intended for Immersion mode.

